



MICRO HAEMATOCRIT CENTRIFUGE

MODEL:KHT-400

SPECIFICATION:

Overall dimension : 25.5 (W) × 23(D) × 24(H) cm
Max speed : 12000 RPM
Max centrifugal force : 14490 × g
No. of tube : 24 places
Timer : 0 - 15 Minutes
Motor : Brush motor
Net Weight : Approx. 7.6 kgs
Power : 110 VAC, 60 HZ or 220VAC, 50HZ (your request)

Note : Capillary tubes of size 75mm in length, 1.5mm O.D.

FEATURES:

- Accurate 15 minutes timer with one minute increment for long runs
- Fast acceleration to 12000 rpm in short time
- Equipped with power-interrupt safety cover
- Comes with safety lock
- 4 large suction feet for installation stability
- Manual electric brake
- Accessory microcapillary tube Reader included.
- The rotor lid for haematocrit testing is easy to install and remove
- Carbon brushes are easily replaced

Ideal for Laboratory of HT, ROUTINE, JAUNDICE in clinic and Veterinary hospital

Specifications subject to change without notice in advance