

# BENCH TOP CENTRIFUGE



## Features :

- The Best Selling and Inexpensive Centrifuge
- Economically priced that everyone can afford
- Durable and salable Centrifuge
- 15 months guaranteed.
- CE marking

## TECHNICAL DATA:

Dimension	270x270x250mm	
Timer	60 minutes	
Choice of motor	Output 90 W	Output 67 W
N.W. (kg)	8.9	7.5
Safety shutout	Option	Option
Standard	CE	CE
Rotor	Various	Various
System	ISO 9002/13488, GMP, CE mark, FDA list	

MODEL	Capacity	Rotor	Max rpm	RCF <sub>xg</sub>	Suggested tube
PLC-03	8x15ml	A-0815	1000 - 4000	2422	Fx12mm conical

## Technical data :

- The **best selling** and less-expensive centrifuge ideal for general purposes in laboratory, hospital, clinic, veterinary, industry and education provided with a desired testing for test tubes containing n vitro diagnostic specimens for qualitative and quantitative test procedures.
- Durable motor incorporating with its carbon brush that can last for 3000 hours in operation.
- Variable speed control from 1000 to 4500 rpm
- Timer activated operation extending to 60 minutes
- Auto balance system to eliminate vibration.
- Balance for rotor and speed setting have been calibrated in plant.
- Choice for additional cover lock fitting and auto shutout safety device for extra costs.
- This equipment complies with EC Electron Magnetic compatibility (EMC) requirement of EN60601-1-2 and IEC801 series regulations.
- This equipment meets the standards of ISO9002, EN46002, ISO13488, CE mark, GMP, FDA.

[\*\*BACK TO PRODUCT INDEX\*\*](#)

[\*\*BACK TO CENTRIFUGE MAIN MENU\*\*](#)