

EXAMINATION LIGHTS

PRODUCT NAME
HALOGEN LAMP

Mobile on castors or wall mountable - supplied with wall plate.

FEATURES

35 watt halogen bulb for whiter, brighter light.

Stainless steel 18/10 central tube.

Reinforced nylon base with 5 feet.

5 Twin wheel swivel castors (50mm).

Flexible 700mm gooseneck arm.

Knock down.

Weighted base.

Power supply: AC Connector - 2 Pin.

CODE

EXA 002 - mobile on castors

EXA 001 - wall type with fastening wall mountable plate

