


**CONTEC™**



**CMS60D-Vet**  
**Pulse Oximeter**





### SpO<sub>2</sub>

- Measurement Range: 0~100%
- Resolution: 1%
- Accuracy: 70~100%, ±2%; 0~69%, unspecified
- Alarm Range: 0~100%

### Pulse Rate

- Measurement Range: 30~250bpm
- Resolution: 1bpm
- Accuracy: ±2 bpm or ±2% (select larger)
- Alarm Range: 0~254bpm

### Storage Environment

- a) Temperature: -40°C ~ +60°C
- b) Relative humidity: 5% ~ 95%
- c) Atmospheric pressure: 500hPa ~ 1060hPa

### Operating Environment

- a) Temperature: 10°C ~ +40°C
- b) Relative Humidity: ≤75%
- c) Atmospheric pressure: 700hPa ~ 1060hPa

### Feature

- Display Mode: 1.8" OLED display
- Power Supply: Two 1.5V (AA size) alkaline battery
- Battery Voltage: DC 2.6V~3.6V
- SpO<sub>2</sub>, Pulse Rate, Bar Graph, Pulse Waveform display
- Battery capacity indication
- Brightness of screen can be changed
- Pulse sound indication

## CMS60D-Vet Pulse Oximeter

- With alarm function
- Synchronous Display on the Oximeter and Computer used USB
- Inner flash memory for 24 hours data storage, and the storage data can be uploaded to computers
- With review function
- Clock function
- It can be connected with an external oximeter probe
- Small in volume, light in weight and convenient in carrying
- **Application** : The device is suitable for being used in family, pet hospital. (It is recommended to use the device when the animal is still.)

**CONTEC™**

CONTEC MEDICAL SYSTEMS CO.,LTD

ADD: #112 Qinhuang West Road, E&T Development Zone,  
Qinhuangdao, Hebei, China  
TEL: +86 335-8015430/8015547  
FAX: +86 335-8015588  
Email: contec88@gmail.com, cms@contecmed.com.cn  
Website: www.contecmed.com