

MRC Handheld Refractometer – REF-312



The REF312 are designed for clinical usage and intended for both veterinary and human beings. Its triple scale reads urine specific gravity, total serum protein and refractive index

REF-312 visual hand held Refractometer with 3 scales and automatic temperature compensation

Scale: Serum protein, urine, refractive index

Measuring range:

- 1) 0 to 12g/dl, resolution 0.2g/dl
- 2) 1.000 to 1.050 SG, resolution 0.002 SG
- 3) 1.3330 to 1.3600 RI, resolution 0.00025 RI

Designed for the measurement of serum protein.

Rugged, durable construction - made from corrosion resistant brass.
High contrast demarcation line for easy, accurate readings.
Easily adjusted zero line by means of readily accessible screw.