

CPWplus L

Large Vet Scale

Durable and portable, CPWplus L series accommodates large animals up to 300kg. The CPWplus L has a removable stainless steel pan that is simple to clean. The backlit display is large and easy to read, and features a hold function to keep results shown on the display.

Model L	CPWplus 150L	CPWplus 200L	CPWplus 300L
Capacity	150kg	200kg	300kg
Readability	0.05kg	0.05kg	0.1kg
Tare Range	Full range		
Pan Size	L - 900x600mm		
Weighing Units	kg, lb, lb:oz, oz		
Calibration	External calibration		
Display	Backlit LCD with 22mm-high digits		
Power Supply	12VDC 500mA adapter / Internal rechargeable battery		
Operating Temperature	0° to 40°C		
Overall Dim.	900x600x80mm (wxdxh)		

#700200059 - Optional rubber mat for CPWplus L

Features

- Vivid, backlit LCD easily visible in any lighting conditions
- Programmable backlight can be set to "always on", "always off" or "light only when weighing"
- Hold function freezes the displayed weight, allowing time to lapse without losing the result
- Selectable digital filtering helps minimise effects of vibration and disturbances
- Rugged construction stands up to industrial use
- Splash-proof housing, keypad and display mean easy cleaning and peace of mind
- Large, grade 304 stainless steel platform allows easy cleaning
- Wall-mount bracket for placing indicator on vertical surface
- Auto power-off to save energy
- Wheels and handle make for easy transportation

