

Specification

Battery life	Lithium Ion 4~4.5 hours continuous use	USB 2.0	For connection to PC & Printer, download images to PC
Viewing options	I-Scan® Video goggles	Video	Composite Video (PAL&NTSC)
Frequency	2.0MHz to 10.0 MHz	Size	40 mm × 250 mm × 120mm
Digital channels	32	Weight	Only 0.6 Kg with battery and TFT LCD.
Cine-Loop	420~450	Temperature range	-10 °C to 40 °C / 14 °F to 104 °F *
Display	5.8 inches Color Transflective TFT	Manufacturer	BMV Technology Co., Ltd

Standard Configuration:

- 1*FarmScan® L60 machine
- 1* Abdominal Convex Probe C350 of 3.5 / 5MHz - (Ewes, goats, sows)
- 1* Li-ion high-power-capacity battery (5 hours/battery)
- 1*Battery charger (AC:110V~250V,50/60Hz)
- 1*Neck strap ,hand strap
- 1*Aluminum alloy carry case
- 1*User guide
- 24 months warranty for machine and 14 months warranty for probes



Accessories – professional your FarmScan® L60



C350– Convex probe



C520– Microconvex probe



LH75– Linear probe



leather protector
(easy-clean by running water)



RC40V– Rectal convex probe



Rectal linear probe



I.R.F™ Introducer



I-Scan® V2
Video Goggles

BMV Technology Co., Ltd

7-8/F Main Building of BMV, ChuangYe 1st Rd,
Qianhai Center, Shenzhen 518100, China
TEL/FAX:(0086)755-26564580
www.bmv.cc info@bmv.cc

*Shaping discovery.
Improving life.*

Specifications subject change without notice

The versatile ultrasound machine for large animal scanning

STANDARD
2
YEAR WARRANTY

Technology Driven
2-Year Warranty



Applications

FarmScan® L60

Thanks to 6 different probes , I.F.R™ introducer and I-Scan® video glasses, FarmScan® L60 can perform all veterinarian applications and even more easily. FarmScan® L60, featuring premium image quality, considerate design for comfortable, fast and reliable detection of pregnancy and gynecological diagnosis in difficult field conditions on a daily basis.

- Fast and reliable detection of pregnancy from 28th day on cows, 13th on mares, 25th on sheep and goats ,18th on pigs (pregnancy control, ovarian examinations, post partum diagnosis, ...)
- Foetal sexing on cows , mares, sheep and goats
- Back fat measurements on pigs, cows and sheep...

Visit www.bmv.cc for more information

[facebook.com/bmvtechnology](https://www.facebook.com/bmvtechnology) | twitter.com/bmvtechnology



Maximum comfort while you work

- Only weight 0.6 kg;
- Excellent imaging quality and ergonomic design;
- Provides versatile solutions with the greatest ease of use for all veterinarian applications – bovine , equine, swine, ovine , Lama and small ruminants;
- 5.8 inch WVGA LCD monitor, wide viewing angle, LCD monitor can turn off when use I-Scan® Video Goggle
- 4.5~5.5 hours Li-ion battery for operating;
- Support USB 2.0 for ultrasound station to transfer images and video to your computer, 450~550-frame cine loop insider;
- Software & Report for reproductive system, and measurement for distance, area, circumference, volume, angel, heart rate;
- Uses a 3rd generation innovative navigation, direct access keys and user buttons for fast, single handed menu selection and setting changes;

The ultrasound scanner that makes you a specialist



I-Scan® video glasses eases your job and gets you focused (optional)



I.F.R™ introducer engineered to lessen fatigue and stress by allowing arm-free scanning, use in a variety of reproductive diagnostic procedures. (optional, compatible with the Rectal convex and Rectal Linear transducers)



Special design rectal convex probe based on customer feedback

Take L60 with you anytime and anywhere



- 10 hours continues working time(with 2 batteries)
- Robust, hard-wired probe designed for everyday use in the difficult field conditions



Leather sheath(easy-clean by running water), designed to be resistant to extremely harsh farm conditions;



Portability by a shoulder belt and hand belt giving you total freedom of movements

Good things come in small packages

- Evaluate the reproductive tract
- Foetal structures
- Foetal sex

Superior Imaging Quality

