



CE 0434



ISO 13485




### Quote from the CEO of SonoScape Co., Ltd.

“ SonoScape is a Chinese leader specializing in developing and manufacturing cutting-edge color Doppler ultrasound systems and ultrasound transducers. Our corporate mission is to deliver revolutionary high-value products and services that improve healthcare worldwide. As the only China Top Brand in the medical industry recognized by General Administration of Quality Supervision, Inspection and Quarantine of the People's Republic of China, SonoScape will be continuing providing the highest quality ultrasound systems and services at the most competitive prices for our customers worldwide.”

*Mr. Jinzhong Yao*

### Quote from the analyst of



“ SonoScape Co. Ltd. has exhibited exceptional entrepreneurial capabilities in developing S8 - the portable colour Doppler System. The advanced 4D capabilities of the ultrasound unit are customised for the high-end European market and yet available at much affordable prices. The company's priority in maintaining high levels of quality in its product lines is well supported by the negligible after-sales costs associated with its ultrasound systems. Also with the launch of SSI-1000, the hand carried ultrasound with multiple clinical applications; SonoScape has made significant strides towards establishing itself as a global medical equipment manufacturer.”

*Mr. Krishanu Bhattacharjee*

# A6V

## Best Choice for Veterinarians



**SonoScape**

Yizhe Building, Yuquan Road, Shenzhen, 518051, China  
Tel: +86-755-26722890 Fax: +86-755-26722850  
E-mail: [sonoscape@sonoscape.net](mailto:sonoscape@sonoscape.net) [www.sonoscape.com](http://www.sonoscape.com)

**SonoScape**

[www.sonoscape.com](http://www.sonoscape.com)

SonoScape is a high-tech company that specializes in manufacturing and developing ultrasound systems and transducers. Supported by numerous innovative technologies, SonoScape introduces A6V as a professional portable ultrasound system for veterinary applications.

- + Elegant outline, only 6 kg weight and optional mobile trolley, A6V creates a comfortable working environment.
- + 12" angle adjustable LCD, two transducer sockets, A6V minimizes your working fatigue.
- + Premium B/W ultrasound technologies, abundant software packages and professional veterinary transducers, A6V assists you in reaching new heights in veterinary applications.



Exclusive 12" adjustable LCD



Thanks to the cutting edge technology, A6V represents our latest designing concept of making easy-to-use ultrasound system in elegant compact design.




- + Graphic Exam Icon: easy to get into exam interface
- + Anti-glare LCD display: displaying uncompromising crystally clear imaging either from distance or under the sunshine
- + One Button quick operation: One-key image optimization, Save, Report, Print, File, Comment, Patient, THI, Clipboard and 9 quick keys for reproductive measurements
- + Abundant veterinary body marks
- + Reliable veterinary design: strengthened case, anti-dust, and anti-splash




To be a HCU dedicated to veterinary applications, A6V makes the ultrasound diagnosis much more reliable and efficient.


- + Veterinary measurement and calculation software packages: early pregnancy diagnosis, cardiac diagnosis, general abdominal diagnosis, reproductive evaluation, meat quality evaluation
- + Extensive veterinary applications: Canine, Hog, Equine, Bovine and Ovine
- + Professional veterinary transducers: anti-corrosion, 3 meters cable, 12 Mhz Linear transducer and 6.5 MHz micro-convex transducer
- + Battery supports for fields usage
- + Super documentation abilities: USB storage, PDF report, AVI/JPEG, DICOM3.0, etc.




Canine Heart



Joint



2B Tendon



Clip Board